

# Technical Data Sheet

## IPTG

for biochemistry

Order number: 1122

Isopropyl- $\beta$ -D-thiogalactopyranoside, or IPTG, is an artificial, non-metabolizable inducer of the *lac* operon<sup>1</sup> of *E. coli*. IPTG causes the *lac* repressor to detach from the operator, allowing transcription of genes in the *lac* operon. These genes include the *lacZ* gene, which encodes  $\beta$ -galactosidase<sup>2</sup> and is probably the most widely used reporter gene in molecular biology research.

Uptake of IPTG into *E. coli* cells occurs via various transport pathways and does not depend (solely) on the  $\beta$ -galactoside permease (*lacY*).

IPTG is soluble in water as well as in organic solvents such as ethanol (20 mg/ml), methanol, DMSO and DMF. For molecular biology applications, the stock solution should be prepared in water to avoid possible physiological effects of the solvents on the cells.

### Application

Since IPTG abolishes the inhibitory effect of the *lac* repressor, it is used for the controlled induction (i.e. selective "switching on") of the transcription of genes that are under the control of the *lac* promoter. Often, these are gene products composed of the actual target gene and  $\beta$ -galactosidase (as a reporter). *E. coli* clones that carry the *lacZ* reporter gene construct (genome- or plasmid-encoded) can be easily selected by blue/white screening (using X-Gal as substrate). The activity of  $\beta$ -galactosidase can be quantified in an enzyme assay using ONPG as a substrate.

Stock solution: 0.1 M (23.8 mg/ml) in water; sterile filtered

Working concentration: 0.1 - 1 mM

### Storage

IPTG is very stable as a solid, even at room temperature. For long-term storage, we nevertheless recommend a temperature of -20°C. IPTG stock solutions should be stored in aliquots at -20°C to avoid repeated freeze-thaw cycles.

JB13072021

<sup>1</sup> The *lac* operon encodes the genes required by *E. coli* for the uptake and degradation of lactose. Naturally, transcription of the *lac* genes is induced by allolactose, a lactose isomer. IPTG, like allolactose, binds to the *lac* repressor, releasing the tetrameric repressor from the *lac* operator. However, unlike allolactose, IPTG cannot be metabolized by the cells, which causes the IPTG concentration in the cell to be kept constant.

<sup>2</sup>  $\beta$ -galactosidase is a hydrolase. It cleaves lactose into galactose and glucose and also converts lactose into allolactose. Since  $\beta$ -galactosidase can also hydrolyse artificial substances, such as ONPG and X-Gal, releasing colored end products, the *lacZ* gene is very often used as a reporter gene in molecular biology.

