



### Challenge

Accurate and reproducible results for real-time PCR experiments.

### Solution

Real-time PCR cyclers with best temperature properties and highest optical precision.

## Highest precision with 1.5-fold discrimination for qPCR quantification applications

### Introduction

A valid real-time PCR assay is characterized by four major properties: PCR efficiency, reproducibility, specificity and sensitivity. Several factors have an influence on these critical points; beside primer design, plastic ware and master mix, of course the performance of the real-time PCR cycler itself affects the results. For this reason, it is mandatory to guarantee best properties of the real-time PCR cycler to rely on the outcome of the applied real-time PCR.

One of the most challenging requirements is to discriminate between small dilution steps, even down to 1:1.5, which proves the dependability and accuracy of the performance of the real-time PCR cycler and therefore of all data created with the device.

## Materials and Methods

For this experiment primers specific for 16S rRNA gene of *E. coli* have been used for amplification of a 120 bp fragment in seven dilution steps with six replicates á 20 µl reaction volume for each concentration. The used master mix was innuMIX qPCR SyGreen Sensitive by Analytik Jena. The time-temperature-profile starts with 2 min at 95 °C, followed by 45 cycles with 5 sec at 95 °C, 5 sec at 58 °C and 15 sec at 72 °C and ends with a melting curve (60 – 95 °C, 15 sec equilibration and 1 °C increment).

## Instrumentation

Data was obtained in color module 1 (470 nm/520 nm, gain = 5) of qTOWER<sup>3</sup> G *touch* for SyGreen.

## Results

The linear regression of the 1.5-fold dilution series results in a PCR efficiency of 97% (slope = -3.41) and  $R^2 = 0.998$ .




Figure 1: Amplification plot for 6-time 1.5-fold dilution of *E. coli* genomic DNA, amplified using primers specific for 16S rRNA gene. With automated threshold settings results are as follows for Ct (SD):  $10^6$  copies, dark red = 16.11 (0.04);  $6,67 \cdot 10^5$  copies, light green = 16.83 (0.05);  $4,44 \cdot 10^5$  copies, dark blue = 17.31 (0.04);  $2,96 \cdot 10^5$  copies, green = 18.05 (0.03);  $1,98 \cdot 10^5$  copies, red = 18.59 (0.03);  $1,32 \cdot 10^5$  copies, yellow = 19.18 (0.03);  $8,78 \cdot 10^4$ , blue = 19.72 (0.03).




Figure 2: Melting peaks of the PCR amplicons result in a melting temperature of 78.80 °C with a SD of 0.13 °C.

## Conclusion

The Analytik Jena qTOWER<sup>3</sup> G *touch* is able to detect and clearly differentiate between minute concentration changes. In the performed 1.5-fold dilution series, high precision discrimination was demonstrated in a reaction with a PCR efficiency of 97% and a  $R^2$  of 0.998.

This document is true and correct at the time of publication; the information within is subject to change. Other documents may supersede this document, including technical modifications and corrections.

### Headquarters

Analytik Jena GmbH  
Konrad-Zuse-Strasse 1  
07745 Jena - Germany

Phone +49 3641 77 70  
Fax +49 3641 77 9279

info@analytik-jena.com  
www.analytik-jena.com

Version: 1.2 | Author: Mja  
en · 12/2020

© Analytik Jena GmbH | © oksix - Fotolia.com