



## Certificate of Quality

[www.qpore.com](http://www.qpore.com)

### Product Information

Part Number: QPPESSFLL033010S  
Lot Number: 10893  
Pore Size: 0.1µm  
Expiry Date: 2021.03.19

### ISO 9001 Quality Standard

The product was manufactured according to ISO 9001 Quality Systems Standard.

### Non-Fiber Releasing

This product meets the criteria for “non-fiber” releasing filters as defined in 21 CFR 210.3 (b) (6)

### Sterilization

This product is gamma irradiate sterilization

### Biosafety

This product meet the criteria of the USP<88>Biological Reactivity test for Class VI Plastics

### Bacterial Endotoxin

Product was determined by LAL test according to current USP. An extract from the product contained was equal to or less than 0.25EU/mL

### Bubble Point

Products were tested according to an established procedure to determine the bubble point of the product. The bubble point of this lot is 26.1-29.0 psi.

### Flow Rate

Products met a flow rate of 34-69ml/min @10psi at 25°C

### Burst Pressure

The burst pressure of products are equal to or more than 87psi.

Tested By:  09-05-2018

Toll free: 1-866-528-4572

Fax: 1-732-412-4040

Website: [www.membrane-solutions.com](http://www.membrane-solutions.com)

Address: JP Morgan International Plaza Suite 650, 14241 Dallas Parkway, Dallas, TX, US

# Certificate of Quality Assurance

Membrane Solutions  
Quality Assurance

---

QC Technician

---

Toll free: 1-866-528-4572

Fax: 1-732-412-4040

Website: [www.membrane-solutions.com](http://www.membrane-solutions.com)

Address: JP Morgan International Plaza Suite 650, 14241 Dallas Parkway, Dallas, TX, US