

PP Reagent Bottles

I. APPLICATIONS

1. Sample storage
2. Dispensing samples

II. FEATURES Material: Strong, medical-grade polypropylene, Autoclavable

1. Made in CleanRoom
2. DNase & RNase free, non-pyrogenic

III. DESIGN

1. High purity
2. Sealing test: Complies with the International Air Transport Association (IATA) Safety Standards
3. Clear bottles enable convenient observation inner content,
4. Ideal choice for storage, transportation and packing of liquid, with uniform thickness of wall, sturdy and durable, which can prevent cracks
5. User-friendly capping/decapping design with proper capping/decapping torques, widely used in industrial automatic filling

IV. DESCRIPTION

Item	Volume	Color	Packaging	Sterile
4-0142	15ml	Clear	50 Stück/Verpackung	No
4-0143	30ml	Clear	50 Stück/Verpackung	No
4-0144	60ml	Clear	25 Stück/Verpackung	No
4-0145	125ml	Clear	10 Stück/Verpackung	No
4-0146	250ml	Clear	5 Stück/Verpackung	No
4-0147	500ml	Clear	5 Stück/Verpackung	No
4-0148	1000ml	Clear	2 Stück/Verpackung	No

V. CERTIFICATION

Moonlab products are manufactured completely according to the TUV ISO 9001:2015 stand

VI. PRODUCT PHOTOS




VII. PRODUCT PARAMETER

Dimensions=mm

Volume	Bottle O.D. D1 (mm)	Bottle height H1 (mm)	Neck I.D. D2 (mm)	Cap O.D. D3 (mm)	Cap height H2 (mm)
15ml	24.75	56.3	14.3	24.1	15
30ml	34	59	21.3	32.2	16.3
60ml	38.7	82.2	21.3	32.2	16.3
125ml	50	95.7	28.3	42.2	17
250ml	59.8	128.3	33.2	46.5	16.9
500ml	72	162.6	43.7	57.3	17.5
1000ml	90	193.9	52.6	68.2	17.5

Volume	Capped bottle height (mm)	Cap spacing (mm)	Max. volume (ml)	Sealing torque (N)
15ml	58.8	1.1	17.5 ml	0.8 N
30ml	62.6	1.15	39 ml	1.0 N
60ml	85.7	1	72 ml	1.0 N
125ml	100	1.75	142 ml	1.2 N
250ml	132.5	2.25	285 ml	2.0 N
500ml	167	1.25	535 ml	2.0 N
1000ml	197.6	0.75	1040 ml	2.5 N



2

