

Specifications

Model	WS-600R
Volume (L)	83
Temperature Range (°C)	15-60
Temperature Accuracy (°C)	±0.1 (38°C)
Temperature Stability (°C)	±1.0 (38°C)
Temperature Controller	Touch Control; Digital Display with PID automatic regulation
Temperature Heater	Heating Plate, 600W
Temperature Timer	Start / Stop / Duration setting; 10sec ~999hr 59min 59sec
Shaker Controller	PID automatic regulation
Shaker Power Unit	Brushless motor
Shaker Methods	Clockwise, Counter-clockwise, forth-back, interrupted,
Shaker Shaking Throw (mm)	20, 30 (standard), 40
Shaker Speed (rpm)	10 ~ 300
Shaker Stability (rpm)	±1 (150rpm)
Platform Dimensions (W × D in mm)	410 x 410
Inside Dimensions (W × D × H in mm)	510 x 510 x 320
Outside Dimensions (W × D × H in mm)	540 x 850 x 625
Net Weight (kg)	102
Max. Load (kg)	15
Order No.	23201