



neoLab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



SERVA SERVAGel™ TG PRiME™ 4 - 12 % ready-to-use gel, 12 sample bags

Price on
request
plus VAT &
Shipping

Product Images



Description

The SERVAGel™ TG PRIME™ 4 - 12 % ready-to-use gel has a proprietary modified Tris/glycine buffer system, which not only guarantees a significantly extended shelf life compared to classical Tris/glycine gels according to Laemmli, but also allows short electrophoresis times when using standard Tris/glycine sample and run buffers (2 gels, 35 min at 350 V const.). It can be used under native and denaturing conditions. The separation range is from 30 up to 300 kDa.

Additional Information

No.	SV-0143
Manufacturer (Brand)	SERVA
EAN	4058072372804
Transport temperature	Wet ice
Number of wells	12
Thickness	1 mm

