



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## SERVA SERVAGel™ TG PRIME™ 4 - 20 % ready-to-use gel, 12 sample bags

Price on  
request

plus VAT &  
Shipping

### Product Images



## Description

---

The SERVAGel™ TG PRIME™ 4 - 20 % ready-to-use gel has a proprietary modified tris/glycine buffer system that not only guarantees a significantly extended shelf life compared to classical Laemmli tris/glycine gels, but also allows short electrophoresis times when using standard tris/glycine sample and run buffers (2 gels, 35 min at 350 V const.). It can be used under native and denaturing conditions. The separation range is from 6 up to 200 kDa.

## Additional Information

---

No.	SV-0147
Manufacturer (Brand)	SERVA
EAN	4058072372842
Transport temperature	Wet ice
Number of wells	12
Thickness	1 mm

