



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



SERVA SERVAGel™ TG PRiME™ 8 % Ready-to-use gel, 10 sample bags

€158.00
plus VAT &
Shipping

Product Images



Description

The SERVAGel™ TG PRIME™ 8 % ready-to-use gel has an in-house developed modified Tris/glycine buffer system, which not only guarantees a significantly extended shelf life compared to classical Tris/glycine gels according to Laemmli, but also allows short electrophoresis times when using standard Tris/glycine sample and run buffers (2 gels, 35 min at 250 V const.). It can be used under native and denaturing conditions. The separation range is from 40 up to 250 kDa.

Additional Information

No.	SV-0124
Manufacturer (Brand)	SERVA
EAN	4058072372613
Transport temperature	Wet ice
Number of wells	10
Thickness	1 mm

