

neoLab Migge GmbH Rischerstr. 7-9 69123 Heidelberg Deutschland +49 (0)6221 / 8442-44 https://www.neolab.de Umsatzsteuer-Identifikationsnummer: DE 143 450 657



## Gilson DISPENSMAN 5 ML (0.5-5ml)

## **Product Images**



## **Description**

DISPENSMAN<sup>™</sup> - precise and safe bottle-top dispenser for many common laboratory fluids DISPENSMAN<sup>™</sup> features an innovative dispensing tap with three position settings. It allows for controlled dispensing of fluids, lossless flushing of the system with fluid, and removal of air in the system. It also has an anti-drip safety feature that returns fluid to the bottle and does not require a protective cap. Volume changes can be adjusted quickly and easily with one hand by pressing to unlock the regulator and sliding it to the desired volume - keeping the bottle securely closed, minimizing the risk of contamination. DISPENSMAN is easily calibrated by the user and is suitable for most liquids, including many acids, bases and solvents. The system is autoclavable and requires no consumables such as pipette tips. Five DISPENSMAN models cover different volume ranges from 0.25 ml to 50 mL. All meet ISO 8655 specifications, and a variety of adapters ensure the DISPENSMAN can be used with most common bottle types.

## **Additional Information**

No.	CF-0033
Manufacturer (Brand)	Gilson
old neoLab article no.	7-4106
EAN	4058072122690
Volume MAX	5 ml
Volume MIN	0.5 ml
Display (reading, setting)	analog
Dispenser properties	mit drehbarem Ventilblock (360°) mit Rückdosierventil mit Adapter
Functionality	mechanisch
Model Bottle Top Dispenser	Dispensman 5 ml
Division	0.1 ml
Material valve spring	Edelstahl
Type Bottle Top Dispenser	Dispenser (analog)
for thread/DM	GL 32 GL 45
adjusted to	"ex"

