



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Socorex Acura® electro 956 Electronic 8- channel microliter pipette, volume 2.5-50 µl

**€1,272.44**  
**plus VAT &**  
**Shipping**

## Product Images

---



## Description

---

**Socorex Acura® electro multichannel pipettes 956 Multi** are flexible electronic pipettes that make pipetting safer than ever.

**Acura® electro Initial Packs** include an electronic pipette, individual control certificate, operating instructions, charging station, batteries with adapter and cable, and samples of Qualitips® pipette tips.

The electronic pipettes feature intuitive software that provides quick access to the six work to the six working functions and the different volume ranges.

- Optimal ergonomics, low weight
- Easy to understand, logical program sequences
- Large display, reversible right / left hand
- Quickly replaceable handle battery
- Fast charging, long working autonomy
- Adjustable ejection for most common tips
- Pipetting cycle counter
- Interchangeable volume units, easy conversion by the user
- The housing of the multichannel pipette can be rotated through 360
- Modular system, convertible with single channel volume unit

Two-year warranty.

## Additional Information

No.	SR-0035
Manufacturer (Brand)	Socorex
old neolab article no.	7-7231
EAN	4058072365868
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	blue
Volume MAX	50 µl
Volume MIN	2.5 µl
sterile	No
Operation via	Knobs/Buttons
Microliter pipettes features	with certification with titration function with mains/charging adapter with charging station with dispensing function including batteries/battery ergonomically shaped IVD compliant
Number of channels	8
Energy supply	Battery
Functionality	electric
Manufacturer	Socorex
Scope of delivery	An electronic pipette, an individual control certificate, operating instructions, a charging station, batteries with adapter and cable, and samples of Qualitips® pipette tips.
Measurement deviation systematic in %	0.8 %
Measurement deviation systematic in µl	0.4 µl
Measurement deviation random in %	0.4 %
Measurement deviation random in µl	0.2 µl
Model Microliter pipette/ Dispenser	Socorex Acura electro 956
Division	0.25 µl
Type Microliter pipette (design)	Air cushion pipette
Type pipette	Microliter pipettes

