




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neoFroxx BSA blocking solution for immunodetection

€188.50 As low
as
€166.95
plus VAT &
Shipping

Product Images

Description

BSA - Blocking Solution for Immunodetection, 500 ml

Western blot membranes must bind proteins and antibodies optimally so that even low protein concentrations can be detected effectively. However, due to the binding property of the membrane surface, insufficient blocking of the free binding sites after protein transfer leads to non-specific background signals and erroneous results. For clean blots and reproducible, reliable results, suitable blocking is therefore mandatory. Depending on the goal of the application (routine; quality control; publication; diagnostics) and the type of sample, different types of blockers are available. For many standard applications, the inexpensive milk powder or albumin (BSA) is the blocker of choice. However, for high background and non-specific bands, a switch to more potent alternatives is recommended.

Additional Information

No.	1550ML500
Manufacturer (Brand)	neoFroxx
EAN	4066306013614
HSNumber	38220000
autoclavable at 121°C	No
Pack	500 ml
Storage	Cooling
sterile	No

