



Thermo Scientific Orion Silver/Sulfide Electrodes

Product Images

Description

Orion Silver/Sulfide Electrodes Thermo Scientific Orion silver/sulfide electrodes can be used to determine silver and sulfide ions in aqueous solution. Connection types: BNC connector Description: solid half cell with BNC connector Deviation: drift-free potentials Diaphragm: free flow diaphragm provides stability Length (metric) Cable: 1m For use with (instruments): ISE meters Operating temperature (metric): 0 ° to 80 °C Type: Ion selective half cell electrode, BNC connector Unit Size: Each

Additional Information

No.	LD-2977
Manufacturer (Brand)	Thermo Scientific

