



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Bandelin® SONOREX ES 4 insert strips

**€99.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

The insert strip is used to hold the insert beakers. When using chemically aggressive liquids (e.g. TICKOPUR J 80 U), removing chemically aggressive or abrasive contaminants or sonicating sample liquids, insert cups are used (indirect sonication). Because nothing may be placed on the bottom of the ultrasonic bath, BANDELIN offers suitable perforated lids to accommodate the insert beakers. Another advantage: When using an insert strip, several insert beakers can be sonicated with different liquids at the same time (e.g. removal of dental cements with STAMMOPUR Z in one and cleaning with simultaneous chemical disinfection of dental instruments with STAMMOPUR DR 8 in a second insert beaker). Please note that in order to transfer the ultrasound to the insert cup, the surface tension of the water must be reduced and therefore a cleaning solution (e.g. TICKOPUR R 33 - 3%) must be used as coupling liquid / contact liquid.

## Additional Information

---

No.	BL-0287
Manufacturer (Brand)	Bandelin
VGKL number	960736010
Transport temperature	Room temperature

