



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Bandelin® SONOREX SH 7 screen holder

€57.00
plus VAT &
Shipping

Product Images



Description

For cleaning test sieves in the ultrasonic bath; stainless steel During ultrasonic operation, nothing must be placed on the bottom of the tank, as abrasion may occur between the sonication material and the ultrasonic oscillating tank. The use of a suitable basket to hold the sonication material is therefore recommended. The BANDELIN baskets are designed in such a way that the sonication material is always in the optimum range of the cavitation field (1-2 cm from the floor) in the oscillating tank. With the aid of a basket, the sonication material can also be conveniently placed in the unit and removed after operation, e.g. for rinsing, without the user coming into contact with the sonication liquid in the oscillating tank.

Additional Information

No.	BL-0121
Manufacturer (Brand)	Bandelin
VGKL number	960732000
Transport temperature	Room temperature

