



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Bandelin® SONOSHAKE SA 1028 shaker attachment

€3,093.00
plus VAT &
Shipping

Product Images



Description

The SA 1028 shaker attachment can be retrofitted to an existing DT 1028 F. The SONOSHAKE offers a wide range of applications for sample preparation in many areas of analysis, e.g. in environmental and food analysis as well as in medical diagnostics. The bath has a footprint of 500 mm x 300 mm, a well depth of only 65 mm and is therefore ideal for sonication of samples in laboratory flasks. The samples can be sonicated in a defined time as well as in continuous operation. Rapid degassing via DEGAS function is also possible. The shaking attachment allows gentle to strong horizontal movement up to a maximum of 20 mm by means of four different shaking frequencies. Both processes can be carried out simultaneously and also separately. For example, prehomogenization can be achieved with a defined shaking frequency and final homogenization with ultrasound in a significantly shorter time. Mains connection: 230 V EU plug

Additional Information

No.	BL-0493
Manufacturer (Brand)	Bandelin
VGKL number	100733249
Transport temperature	Room temperature

