



DURAN[®] Viton gasket, 90 mm

Product Images



Description

FKM seals, Ø 90 mm One of the most important separation methods in the laboratory environment is filtration and is used for substance separation or purification of different media. Filters are required for this working method, which are selected according to the application in terms of size and porosity. The filter attachment consisting of a receiver vessel, sealing ring, glass filter plate and funnel can be disassembled into its individual parts and each individual component can be replaced or varied to adapt the filter to the application. Thus, the three unit sizes combined with four filter plates each of different porosity create 12 different filter options. Special features and advantages: Due to the seal, a tight, leak-proof connection of the individual components of the filter attachment is possible Burned-in and therefore permanently durable white printed image Application: Gasket

Additional Information

No.	SD-2702
Manufacturer (Brand)	DURAN
VGKL number	112922090

