



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Fritsch® impact bar P-14 for pre-cutting/difficult specimens

€696.00
plus VAT &
Shipping

Product Images

Description

By using the beater bar, a sieve ring with a larger diameter and a rotor, the material to be ground is comminuted with increased shearing action. This enables easy pre-crushing or fine crushing of hard-brittle to soft samples, temperature-sensitive and fat-containing samples, as well as samples with residual moisture. For the discharge of large sample quantities, the conversion kit for comminution of larger quantities is recommended. Please note: For this impact bar you need an impact rotor and at least one special sieve ring for the stainless steel impact bar.

Additional Information

No.	FR-0547
Manufacturer (Brand)	Fritsch
Transport temperature	Room temperature

