



Sartorius® Microsart ATMP Extraction

Product Images

Description

To achieve the highest possible sensitivity and to avoid inhibitory effects during PCR runs, it is recommended to perform a DNA extraction beforehand. DNA extraction kits are available for almost all sample types. However, almost all offered extraction kits are not free of contamination with microbial DNA, as it can be found ubiquitously. Microsart® ATMP Extraction is a unique extraction kit that enables contamination-free DNA extraction from cell cultures and ATMPs.

Additional Information

No.	ST-6544
Manufacturer (Brand)	Sartorius

