



Sartorius® Vivapure LentiSELECT 1000

Product Images



Description

For rapid purification of up to 2.5×10^9 infectious particles per ml from a 500 ml cell culture supernatant. This kit contains all reagents and units necessary for purification of 2.5×10^9 infectious particles from 500 ml cell culture volume. The entire purification procedure is simple and pump-driven, which reduces labor. Unlike the conventional purification method such as ultracentrifugation, Vivapure LentiSELECT 500 provides a fast and easy method for purification of VSV-G pseudotype lentivirus, eliminating the need for expensive purification equipment such as an ultracentrifuge. Pump operated kit Easy and controlled purification of lentivirus from large cell culture volumes. Fast Non-toxic Convenient Cost-effective Better purity than pelleting by ultracentrifugation

Additional Information

No.	ST-6576
Manufacturer (Brand)	Sartorius

