



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Usbeck® Bunsen burner propane, DIN, gas tap,
pilot flame, air regulation**

**€47.70
plus VAT &
Shipping**



Description

according to DIN 30665. DIN burners cannot be shut off completely according to standard 30665. Therefore, a small flame continues to burn even when the gas tap or needle valve is closed. The gas supply must be shut off at the central shut-off in the laboratory

Additional Information

No.	UB-0552
Manufacturer (Brand)	Usbeck
Transport temperature	Room temperature

