



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Rinsing device according to Swiegot for NMR tubes

€177.00
plus VAT &
Shipping

Product Images



Description

NMR tubes must always exhibit a high degree of cleanliness and purity. To ensure that the NMR measurement results are not falsified, the tubes should be rinsed thoroughly beforehand.

The neoLab rinsing **device for NMR tubes** offers the optimal solution for these high demands.

- no more time-consuming and tedious manual rinsing necessary
- lower consumption of solvents for the rinsing process
- Attachment (100 ml) with 500 ml bottle, with PE olive. Suitable for vessels up to 20 mm diameter.

Additional Information

No.	E-3993
Manufacturer (Brand)	neoLab
EAN	4058072161194
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Material	Borosilicate glass
Form	straight
Connection General / Tempering jacket / Cooler / Vacuum	Olive
Standard ground core	29/32
Type connection	Hose olive
Type glass jar	Rinsing device for tubes
Type ground	Standard ground

