



## Nalgene® Mikroröhrchen mit Schraubverschluss

### Product Images

### Beschreibung

Mikroröhrchen mit Schraubverschluss Echter Schutz : sichere Vorbereitung, Übertragung und Lagerung von Proben mit Thermo Scientific Mikroröhrchen mit Schraubverschluss. Erhältlich in mehreren Volumen, Verpackungsoptionen und Reinigungsangaben. Beschreibung: 0.5ml Screw Cap Tube, Knurl, Skirted, Natural, Non-Sterile tube w/ unattached cap Autoklavierbar: No Kapazität (metrisch): 0.5 mL Farbe: Naturfarben Max. Kraft: 18000xG Verpackung: 500 x pk, 10 pk x cs, 5000 x cs. Single Bagged. ZipLock bag. Temperaturbereich (metrisch): Cold Storage: -80deg; to plus; 95deg; C. Liquid Nitrogen: Short term Vapor Phase Quick Flash Freeze acceptable, not recommended for LN2 submersion storage Reinigungsnachweise: Non-Sterile Free of: Rnase/Dnase, Human and Mouse DNA, Pyrogen Endotoxin Bindungseigenschaft: Low Nucleic Acid and Protein Binding Low Enzyme Inhibition Max. Drehmoment: Manual/Hand Torque: Standard: 0.125 lbf.in to 2.5 lbf.in; Metric: 1.412 Ncm to 28.24 Ncm. Automated/Robotic Torque Standard: 1.312 lbf.in to 2.5 lbf.in; Metric: 14.83 Ncm to 28.24 Ncm Max. RZB: 18,000 x g Graduated: Yes Unit Size: Case of 5000

### Zusätzliche Informationen

|                    |         |
|--------------------|---------|
| Artikelnummer      | CD-3184 |
| Hersteller (Marke) | Nalgene |

