



Nalgene® Mikroröhrchen mit Schraubverschluss

Product Images

Beschreibung

Mikroröhrchen mit Schraubverschluss Echter Schutz : sichere Vorbereitung, Übertragung und Lagerung von Proben mit Thermo Scientific Mikroröhrchen mit Schraubverschluss. Erhältlich in mehreren Volumen, Verpackungsoptionen und Reinigungsangaben. Beschreibung: 0.5ml Screw Cap Tube, NonKnurl,Skirted,Natural, E-Beam Sterile tube w/ attached cap Autoklavierbar: No Kapazität (metrisch): 0,5 ml Max. Kraft: 25000xG Verpackung: 500 x pk, 10 pk x cs, 5000 x cs. Single Bagged. ZipLock bag. Temperaturbereich (metrisch): Cold Storage: -80deg; to +95deg;C. Liquid Nitrogen: Short term Vapor Phase Quick Flash Freeze acceptable, not recommended for LN2 submersion storage Reinigungsnachweise: Sterile Electron Beam, SAL 10sup-6/sup Free of: Rnase/Dnase, Human and Mouse DNA, Pyrogen, Endotoxin, Bioburden, ATP PCR Inhibitor Bindungseigenschaft: Low Nucleic Acid and Protein Binding Low Enzyme Inhibition Max. Drehmoment: Manual/Hand Torque: Standard: 0.125 lbf.in to 2.5 lbf.in; Metric: 1.412 Ncm to 28.24 Ncm. Automated/Robotic Torque Standard: 1.312 lbf.in to 2.5 lbf.in; Metric: 14.83 Ncm to 28.24 Ncm Max. RZB: 18.000 x g Graduated: Ja Unit Size: Each

Zusätzliche Informationen

Artikelnummer	CD-3076
Hersteller (Marke)	Nalgene

