



## HiMedia 1,4-Dibromobutane, 250 ml

### Artikelbilder

---



# Beschreibung

Hi-Media 1,4-Dibromobutane, 250 ml

## Zusätzliche Informationen

|                    |                |
|--------------------|----------------|
| Art.-Nr.           | RM9431-250ML   |
| ClassPG            | 6.1            |
| Hersteller (Marke) | Hi-Media       |
| HSNummer           | 29033919       |
| Packung            | 250 ml         |
| Storage            | Raumtemperatur |
| UNNummer           | UN2810         |

