



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



HiMedia Fixative (Buffered Formalin fixative) for fixing cytological or histolo

56,90 €
zzgl. MwSt &
Versand

Product Images



Beschreibung

Hi-Media Fixative (Buffered Formalin fixative) for fixing cytological or histological samples, 500 ml

Zusätzliche Informationen

Art.-Nr.	S109-500ML
Hersteller (Marke)	Hi-Media
HSNummer	38220000
Transporttemperatur	Raumtemperatur
autoklavierbar bei 121°C	Nein
Packung	500 ml
Storage	Raumtemperatur
steril	Nein

