



# HiMedia L-Glutamine, Certified, 500 g

## Artikelbilder

## Beschreibung

Hi-Media L-Glutamine, Certified, 500 g

## Zusätzliche Informationen

|                    |                |
|--------------------|----------------|
| Art.-Nr.           | RM049-500G     |
| Hersteller (Marke) | Hi-Media       |
| CASNummer          | 56-85-9        |
| HSNummer           | 29241900       |
| Packung            | 500 g          |
| Storage            | Raumtemperatur |
| steril             | Nein           |

