



## HiMedia L-Methionine, 100 g

### Artikelbilder

---



# Beschreibung

Hi-Media L-Methionine, 100 g

## Zusätzliche Informationen

|                    |                |
|--------------------|----------------|
| Art.-Nr.           | PCT0315-100G   |
| Hersteller (Marke) | Hi-Media       |
| CASNummer          | 63-68-3        |
| HSNummer           | 29304010       |
| Packung            | 100 g          |
| Storage            | Raumtemperatur |
| steril             | Nein           |

