



HiMedia N6-2'-O-Dibutyryladenosine-3',5' cyclic, 25 mg

Artikelbilder

Beschreibung

Hi-Media N6-2'-O-Dibutyryladenosine-3',5' cyclic, 25 mg

Zusätzliche Informationen

Art.-Nr.	TC467-25MG
Hersteller (Marke)	Hi-Media
CASNummer	16980-89-5
Transporttemperatur	Trockeneis
Packung	25 mg
Storage	Raumtemperatur
steril	Nein

