



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



HiMedia Sodium tellurate dihydrate suitable as bacteriological reagent, 100 g

384,80 €
zzgl. MwSt &
Versand

Product Images



Beschreibung

Hi-Media Sodium tellurate dihydrate suitable as bacteriological reagent, 100 g

Zusätzliche Informationen

Art.-Nr.	RM166-100G
ClassPG	6.1
Hersteller (Marke)	Hi-Media
HSNummer	28429090
Transporttemperatur	Raumtemperatur
autoklavierbar bei 121°C	Nein
Packung	100 g
Storage	Raumtemperatur
UNNummer	UN3288
steril	Nein

