




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Whatman® 593 Filtertüten 185 mm 100/Pack

65,53 €
zzgl. MwSt &
Versand

Product Images

Beschreibung

593 FF 185 mm 100/Pk A standard grade filter papers for fine precipitates. Increased filter surface area increases total loading capacity and decreases overall filtration time. Reduced contact between filter and funnel maintains flow rate. Eliminated need to quadrantfold filter circles to fit conical funnels saves time in repetitive or multiple analyses. Prepleated filters have the same technical properties as the respective filter circles. Cellulose filters are used in qualitative analytical techniques to determine and identify materials. Prepleated qualitative filters give improved flow rate and increased loading capacity compared to equivalent quadrant folded filters.

Zusätzliche Informationen

Art.-Nr.	VL-0137
Hersteller (Marke)	Whatman
VGKL Nummer	250206306
Transporttemperatur	Raumtemperatur

