




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Whatman® Sorte 0858 ½ qualitatives Filterpapier, 240 mm

68,57 €
zzgl. MwSt &
Versand

Product Images

Beschreibung

Grade 113V Qualitative Filter Paper Folded (Prepleated), 240 mm Grade 0858 ½ qualitative filter papers Savings in time required to quadrantfold circles to fit conical filter funnels in repetitive or multiple analyses. Decreased overall filtration time because of the extra surface area exposed; the normal slow down of filtration speed due to the loading of particulate is postponed. Increased total loading capacity as more filter area is available. Maintained flow rate due to the reduction in filter paper contact with funnel side and the selfsupporting shape of the filter itself. The prepleating does not significantly affect any of the technical data, and the same figures may be used for the flat circles. Timesaving Whatman qualitative grades are offered in this convenient form, which gives major advantages over flat circles. Grained, with medium fast flow rate and medium retention. A universal filter paper, used for the filtration of extracts, oils, beer, syrups etc.; also applied in filter presses or for the aspiration of liquids. Available flat stock form as Grade 0858. Grade 0858 ½: 712 µm

Zusätzliche Informationen

Art.-Nr.	VL-0248
Hersteller (Marke)	Whatman
Transporttemperatur	Raumtemperatur

