



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



---

## Whatman® Sorte 0858 ½ qualitatives Filterpapier, 320 mm

**108,31 €**  
**zzgl. MwSt &**  
**Versand**

### Product Images

---

# Beschreibung

Grade 113V Qualitative Filter Paper Folded (Prepleated), 320 mm Grade 0858 ½ qualitative filter papers Savings in time required to quadrantfold circles to fit conical filter funnels in repetitive or multiple analyses. Decreased overall filtration time because of the extra surface area exposed; the normal slow down of filtration speed due to the loading of particulate is postponed. Increased total loading capacity as more filter area is available. Maintained flow rate due to the reduction in filter paper contact with funnel side and the selfsupporting shape of the filter itself. The prepleating does not significantly affect any of the technical data, and the same figures may be used for the flat circles. Timesaving Whatman qualitative grades are offered in this convenient form, which gives major advantages over flat circles. Grained, with medium fast flow rate and medium retention. A universal filter paper, used for the filtration of extracts, oils, beer, syrups etc.; also applied in filter presses or for the aspiration of liquids. Available flat stock form as Grade 0858. Grade 0858 ½: 712 µm

# Zusätzliche Informationen

Art.-Nr.	VL-0250
Hersteller (Marke)	Whatman
Transporttemperatur	Raumtemperatur

