




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Whatman® GD/XP 25/0.45 PTFE 150/Pack

554,00 €
zzgl. MwSt &
Versand

Product Images

Beschreibung

GD/XP Syringe Filter, 0.45 µm, PTFE filtration medium Whatman GD/XP Disposable Syringe Filters are highly suitable for use with samples that require inorganic ion analysis as levels of ion extractables are minimized. They are also an alternative choice for users requiring a filter that exhibits extremely low protein binding characteristics. GD/XP Syringe Filters contain a twolayer prefilter stack comprised of 20 µm and 5 µm polypropylene filters. The last stage of filtration is a choice of membrane, which is positioned below the prefilter stack.

Zusätzliche Informationen

Art.-Nr.	VL-3169
Hersteller (Marke)	Whatman
Transporttemperatur	Raumtemperatur

