




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Whatman® IC 25/0.2 200/Pack

3.492,00 €
zzgl. MwSt &
Versand

Product Images

Beschreibung

original description: Anotop 25 Syringe Filter IC, 0,2µm 200/pk (WEB) Anotop 25 Syringe Filter IC, 0.2 µm Anotop 25 syringe filter IC Each batch certified for ion chromatography (IC). Enhanced consistency of analytical results. Extended column life. Certified and guaranteed low levels of anion leaching for improved results. Anotop Syringe Filters are a universal solution for numerous filtration applications. Anotop filters can be used with most organic solvents and aqueous materials and are suitable for sample volumes up to 100 ml. The distinctive hexagonal housing is manufactured from pigmentfree polypropylene to eliminate sample contamination. No wetting agents or adhesives are used in the manufacturing process. Anotop Syringe Filters contain the proprietary aluminabased Anopore membrane and are supplied in three pore sizes. Glass microfiber prefilter versions are available for difficulttofilter samples. Whatman Anotop IC Syringe Filters are specifically designed for the preparation of samples for subsequent ion chromatography and HPLC analysis. These devices ensure very low levels of anion leaching for ion chromatography testing. Anotop IC

Zusätzliche Informationen

Art.-Nr.	VL-3083
Hersteller (Marke)	Whatman
Transporttemperatur	Raumtemperatur

