



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Whatman® Polytetrafluorethylen Spritzenvorsatzfilter, 1 µm, unsteril (500 Stk.)

1.694,00 €
zzgl. MwSt &
Versand

Product Images

Beschreibung

original description: ReZist Syringe Filter 30mm, PTFE, 1,0µm 500/pk (WEB) ReZist Syringe Filter, 30 mm, PTFE, 1.0 µm ReZist syringe filters Hydrophobic PTFE membrane is laminated with polypropylene. Excellent chemical resistance against standard organic HPLC solvents. Permits optimal utilization of small sample volumes. Integral, permanently hydrophobic PTFE membranes. Polypropylene support. Extremely high chemical resistance. Whatman ReZist range of syringe filters has been specifically designed to be resistant to organic solvents. These filters are highly suitable for the clarification of aggressive organic solvents. ReZist 30 mm filters can also be used as a venting filter for small vessels.

Zusätzliche Informationen

Art.-Nr.	VL-0764
Hersteller (Marke)	Whatman
Transporttemperatur	Raumtemperatur

