



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



---

## Whatman® Polytetrafluorethylen Spritzenvorsatzfilter, 5 µm, unsteril (500 Stk.)

**1.694,00 €**  
**zzgl. MwSt &**  
**Versand**

### Product Images

---

## Beschreibung

---

original description: ReZist Syringe Filter 30mm, PTFE, 5,0µm 500/pk (WEB) ReZist Syringe Filter, 30 mm, PTFE, 5.0 µm ReZist syringe filters Hydrophobic PTFE membrane is laminated with polypropylene. Excellent chemical resistance against standard organic HPLC solvents. Permits optimal utilization of small sample volumes. Integral, permanently hydrophobic PTFE membranes. Polypropylene support. Extremely high chemical resistance. Whatman ReZist range of syringe filters has been specifically designed to be resistant to organic solvents. These filters are highly suitable for the clarification of aggressive organic solvents. ReZist 30 mm filters can also be used as a venting filter for small vessels.

## Zusätzliche Informationen

---

Art.-Nr.	VL-0766
Hersteller (Marke)	Whatman
Transporttemperatur	Raumtemperatur

