



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Whatman® MV050 Serie MV 050A/0
Vakuumfiltrationsapparat mit
Schnellverschlussklemme, Edelstahl, 500 ml,
47/50 mm**

**1.298,00 €
zzgl. MwSt &
Versand**

Product Images

Beschreibung

MV050 Series MV 050A/0 Vacuum Filtration Apparatus with rapid closure clamp, stainless steel, 500 ml, 47/50 mm All MV series vacuum filtration devices are made of stainless steel, which is especially suitable for microbiological applications. The system can be used up to 200°C, is autoclavable, and can be sterilized by dry heat up to 180°C.

Zusätzliche Informationen

| | |
|---------------------|----------------|
| Art.-Nr. | VL-0618 |
| Hersteller (Marke) | Whatman |
| VGKL Nummer | 540240020 |
| Transporttemperatur | Raumtemperatur |

