



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



---

## Whatman® Membrane Holder, Glass Support, 90mm, 1000ml volume 1/Pack

**1.333,00 €**  
**zzgl. MwSt &**  
**Versand**

### Product Images

---

## Beschreibung

---

VacuumType Glass Membrane Holder, glass support, 90 mm, 1000 ml volume, membrane filter accessories Membrane holder produced from borosilicate glass, suitable for aqueous and organic solvent filtration. Produced from borosilicate glass and available with a choice of support screen. Suitable for aqueous and organic solvent filtration. The funnel seal ensures that the sample does not bypass the membrane and that particulates are retained on the surface of the membrane. Vacuum type glass holders The sintered glass support is recommended for filtration and biological analysis. The 304 stainless steel support screen is suitable for use with proteinaceous solutions.

## Zusätzliche Informationen

---

Art.-Nr.	VL-2643
Hersteller (Marke)	Whatman
Transporttemperatur	Raumtemperatur

