



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Whatman® Membrane Holder, Stainless Steel Syringe Filter Type, 25mm 1/Pack

401,00 €
zzgl. MwSt &
Versand

Product Images

Beschreibung

Syringe Type Membrane Holder, stainless steel, 25 mm Membrane holder, stainless steel syringe filter type Available in stainless steel and polypropylene with luer fittings for use with a standard syringe. The holders are designed for the quick and easy clarification, sterilization, and removal of particulates from small volume samples, typically for HPLC applications. The holders contain PTFE gaskets and Orings and allow the membrane to be autoclaved in place without the filter sticking to the holder. Luer lock fittings connect to a standard syringe and offer convenience and ease of use for clarification, sterilization, and removal of particulates from small volumes of liquid (e.g., HPLC samples and solvents).

Zusätzliche Informationen

Art.-Nr.	VL-2649
Hersteller (Marke)	Whatman
Transporttemperatur	Raumtemperatur

