




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Thermo Scientific Matrix ScrewTop Tri-coded Tubes

486,00 €
zzgl. MwSt &
Versand

Product Images

Beschreibung

Matrix ScrewTop Tri-coded Tubes Securely store precious samples at low temperatures, including vapor phase Liquid Nitrogen with Thermo Scientific Matrix 500l and 1.0mL ScrewTop Tubes. Kapazität (angloamerikanisch): 0,0338 fl oz. DNase-frei: Ja Typ: Röhrchen mit Schraubgewinde oben Beschreibung: V-Bottom ScrewTop Tubes with white patch and linear barcode Verschlussfarbe: Colorless Kapazität (metrisch): 1.0 mL Bodenform: V-Boden Packungskonfiguration: Bulk Wellform: V-Boden Barcode: Non-proprietary, 12x12 Data-matrix with ECC200 Built-in Error Correction Unit Size: Each

Zusätzliche Informationen

Art.-Nr.	TA-0044
Hersteller (Marke)	Thermo Scientific
Transporttemperatur	Raumtemperatur

