



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Corning® 2 mL Assorted Colors Internal Threaded Polypropylene Cryogenic Vial, Self- Standing with Round Bottom

397,31 €
zzgl. MwSt &
Versand

Product Images



Beschreibung

VIAL,CRYOGENIC,INT-THREAD 2.0mL,RND,SELF STDG,W/WSHR,50/500 (Rainbow cap) Case Size: 10 Pack size: 50

Zusätzliche Informationen

Art.-Nr.	XC-0850
Hersteller (Marke)	Corning
Transporttemperatur	Raumtemperatur

