



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Corning® 2 mL White Cap Internal Threaded Polypropylene Cryogenic Vial, Self-Standing with Round Bottom

362,84 €
zzgl. MwSt &
Versand

Product Images



Beschreibung

VIAL,CRYOGENIC,INT-THRD,2.0mL,RD,SLF-STDNG,W/WSHR,WHT CAPS Case Size: 10 Pack size: 50

Zusätzliche Informationen

Art.-Nr.	XC-0855
Hersteller (Marke)	Corning
Transporttemperatur	Raumtemperatur

