



Macherey & Nagel® NucleoTraPCR (100)

Artikelbilder



Beschreibung

NucleoTraP CR (100) 100 preps for PCR clean-up - NucleoTraP CR Suspension, buffers pack of 100 __UN 3211 Perchlorate, anorganisch, wässrige Lösung, n.a.g. 5.1 III 0,100 L ADR/GGVSEB 01 (E) LQ PAX 551 2.5 L CAO 555 30 L Q: 0.040 / 0.004

Zusätzliche Informationen

Art.-Nr.	ML-4195
Hersteller (Marke)	Macherey & Nagel

