



## Sartorius® Sartolab MultiStation

### Artikelbilder

### Beschreibung

For hands-free parallel filtration of up to 6 samples: Sartolab® MultiStation Sartolab® MultiStation is a stand specially designed to hold 1 to 6 vacuum filtration units, allowing simultaneous filtration of up to 6 samples. The MultiStation is permanently connected to your vacuum source. Easily install your vacuum filtration units in the MultiStation for quick and easy filtration of samples without the need for installation of extra connectors and time-consuming stabilization. Sartolab® MultiStation works with our Sartolab® RF 50 vacuum filtration units; the funnel adapter of these units is designed to fit perfectly in the bracket of the MultiStation. With one click, connect the filtration unit to the device, assuring perfect filter stability. With a second click, engage the vacuum automatically and begin filtering. Easily manipulate your samples with the rotating, multi-directional head, and easily keep track of your samples during filtration with the numbered brackets. User Benefits One single vacuum source enables simultaneous filtration of up to 6 samples Time-saving (no installation time for each filter unit before use) Hands-free filtration

### Zusätzliche Informationen

Art.-Nr.	ST-6765
Hersteller (Marke)	Sartorius

