



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Thermo Scientific Smart-Vue Client Server Software Smart-Vue Software 2.0, 20-99

3.020,00 €
zzgl. MwSt &
Versand

Product Images



Beschreibung

Smart-Vue Client Server Software Automate monitoring of vital laboratory equipment parameters with Thermo Scientific™ Smart-Vue™ Client Server Software. Find peace of mind knowing your work is fully monitored 24/7 with immediate notifications reaching you even when you are away from the lab. This software provides a plug and play solution with an intuitive interface. Intelligent Plug Play Simplicity Installation Wizard for simplified, step-by-step process Service Discovery Protocol (SDP) enables every module to send signal to receiver—identifies most efficient communication route based on RF signal strength Module automatically establishes new path to receiver if no communication during 6-hour period Back-up system with memory capacity up to 3000 points, with uninterruptible power for continuous data logging Intuitive Dashboard Daily usage is simplified Single-click export of reports to Microsoft™ Excel™, Word™ or Adobe™ PDF formats Flexible creation of graphs summarizing data from individual or multiple sensors Visual depiction of RF signal strength an QoS as monitored by each module Data Traceability Security Secure, customized levels of user access Audit trails Archived alarms Password protection Secure digital transmission and data back-ups Numerous Configuration Capabilities Combines RF and TCP/IP network technologies Share equipment data over entire network instantly Download simple plug-in to visualize equipment data readings from any workstation Compatible with: Refrigerators, Freezers, LN 2 Storage, CO 2 Incubators, Clean Rooms, Storage Areas, Ovens Recommended for: Ideal for use in highly controlled environments, including pharmaceutical production, biotechnology, IVF clinics, clinical studies and diagnostic laboratories. Bezeichnung: Smart-Vue Software 2.0, 20-99 sensors Enthält: Software drivers for all system receivers and repeaters; includes user manuals / quick installation guides. Typ: Sensor Typ: Software Zertifizierungen/Konformität: Smart-Vue facilitates conformance to 21 CFR Part 11, cGxP and SOP compliance.

Zusätzliche Informationen

Art.-Nr.	LD-5250
Hersteller (Marke)	Thermo Scientific

