



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Wiggins® PT100 Temperature sensor, Type I Length: 170 mm, Diameter: 4 mm, Materi

239,06 €
zzgl. MwSt &
Versand

Product Images

Beschreibung

Wiggins PT100 Temperature sensor, Type I Length: 170 mm; Diameter: 4 mm; Material: Stainless steel; Admissible temperature: -30 ~ +300°C

Zusätzliche Informationen

Art.-Nr.	WI-0032
Hersteller (Marke)	Wiggins
Transporttemperatur	Raumtemperatur

