



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Thermo Scientific Fiberlite 1000mL Bottles O- ring for 1000mL bottles Pack of 2

16,10 €
zzgl. MwSt &
Versand



Beschreibung

Fiberlite 1000mL Bottles Thermo Scientific™ 1000mL Fiberlite Bottles combine a Nalgene bottle body with newly designed Fiberlite Cap Closure. These bottles are available for Thermo Scientific™ Sorvall™ LYNX Superspeed and Bios A Large-capacity centrifuges Features • PPCO PC • Increased cycle life compared to previous Thermo Scientific bottles - 100 cycle life (caps + bottles) • Certifications USP Class VI - Non-Cytotoxic • Max speeds match rotors • Ergonomic speed handle design allows much easier bottle handling and removal
Bezeichnung: O-ring for 1000mL bottles
Material: Silicon
Verpackung: Set of 4 Volumen bis Rand (metrisch): 1000mL

Zusätzliche Informationen

Art.-Nr.	LD-0001
Hersteller (Marke)	Thermo Scientific
Transporttemperatur	Raumtemperatur

