



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Thermo Scientific Fiberlite 250mL Bottle Polycarbonate (PC) Pack of 2 Set of 6

282,00 €
zzgl. MwSt &
Versand

Product Images



Beschreibung

Fiberlite 250mL Bottle Thermo Scientific™ 250mL Fiberlite Bottles combine a Nalgene bottle body with newly designed Fiberlite Cap Closure. These bottle will be available for Thermo Scientific™ Sorvall™ LYNX Superspeed and General Purpose Centrifuges. Features • PPCO PC • Increased cycle life compared to previous Thermo Scientific bottles - 100 cycle life (caps + bottles) • Certifications USP Class VI - Non-Cytotoxic • Max speeds match rotors • Ergonomic speed handle design allows much easier bottle handling and removal
Bezeichnung: Polycarbonate (PC) Verpackung: Set of 6
Volumen bis Rand (metrisch): 250mL

Zusätzliche Informationen

Art.-Nr.	LD-0009
Hersteller (Marke)	Thermo Scientific
Transporttemperatur	Raumtemperatur

