



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Thermo Scientific Fiberlite 250mL Bottle Polypropylene (PPCO) Pack of 2 Set of 6

**261,00 €**  
**zzgl. MwSt &**  
**Versand**

## Product Images

---



## Beschreibung

---

Fiberlite 250mL Bottle Thermo Scientific™ 250mL Fiberlite Bottles combine a Nalgene bottle body with newly designed Fiberlite Cap Closure. These bottle will be available for Thermo Scientific™ Sorvall™ LYNX Superspeed and General Purpose Centrifuges. Features • PPCO PC • Increased cycle life compared to previous Thermo Scientific bottles - 100 cycle life (caps + bottles) • Certifications USP Class VI - Non-Cytotoxic • Max speeds match rotors • Ergonomic speed handle design allows much easier bottle handling and removal  
Bezeichnung: Polypropylene (PPCO) Verpackung: Set of 6 Volumen bis Rand (metrisch): 250mL

## Zusätzliche Informationen

---

Art.-Nr.	LD-0008
Hersteller (Marke)	Thermo Scientific
Transporttemperatur	Raumtemperatur

