




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Thermo Scientific Owl D2 Wide-Gel- Electrophoresesystem Rapid load 25 (2X) Each 1.5mm

59,20 €
zzgl. MwSt &
Versand



Beschreibung

Owl D2 Wide-Gel-Electrophoresesystem Owl D2 Wide-Gel Electrophoresis System Standard, 40 teeth, 1.0mm thick tooth Eliminate comb assembly with Thermo Scientific™ Owl™ D2 Wide-Gel Electrophoresis System Combs. These heavy-duty, one-piece combs are for use with Owl™ D2 Wide-Gel Electrophoresis Systems. Bezeichnung: Rapid load Zur Verwendung mit: Thermo Scientific™ Owl D2 Wide-Gel Electrophoresis Systems Artikelbeschreibung: Rapid load; 25 (2X) teeth; 1.5mm thick Anzahl der Zähne: 25 (2X) Dicke der Zähne (metrisch): 1.5mm Zur Verwendung mit (Geräte): Owltrade; D2 Wide-Gel Electrophoresis Systems Typ: Combs

Zusätzliche Informationen

Art.-Nr.	LD-3712
Hersteller (Marke)	Thermo Scientific

